

Summary of Substances & Concentrations without Interference

Substance	Concentration Tested (mg/dL)	Therapeutic/ Physiologic Concentration Range (or upper limit) (mg/dL)	FORA GD20
Acetylsalicylic Acid	50	2 - 10	No interference
Ascorbic Acid	4	2	No interference
Acyclovir	3.1	0.23 - 0.31	No interference
Allopurinol	5	0.5	No interference
Amitriptylline	0.25	0.012 - 0.025	No interference
Amoxicillin	11	0.55 - 1.1	No interference
Ampicillin	5	0.5	No interference
Aspirin (Salicylic Acid)	60	10 - 30	No interference
Atenolol	10	0.1 - 0.2	No interference
Bicarbonate	336 (40 mM)	244 (29 mM)	No interference
Bile Acids (Cholic Acid)	6	0.7	No interference
Caffeine	10	0.3 - 1.5	No interference
Calcium	5 mM	2.8 mM	No interference
Chloride	140 mM	108 mM	No interference
Cholesterol	500	300	No interference
Clonidine	2	0.0001 - 0.0002	No interference
Creatinine	5	1.7	No interference
Digoxin	0.16	0.0001 - 0.00025	No interference
Diphenhydramine	1	0.01 - 0.1	No interference
EDTA (K3)	702	175.5	No interference
Enalapril	0.15	0.012 - 0.015	No interference
Erythromycin	20	0.2 - 2.0	No interference
Estrone	0.1	0.0011	No interference
Famotidine	0.13	0.008 - 0.013	No interference
Fluoxetine	0.8	0.08	No interference
Fructose	1000	7.5	No interference
Furosemide	2	0.1 - 0.3	No interference
Galactose	1000	7.5	No interference
Gentisic Acid	2	0.2 - 0.6	No interference
Glyburide	1.07	0.018 - 0.025	No interference
Hemoglobin	500	2.5	No interference
Heparin (Na)	6800 U/dL	4000 U/dL	No interference
Ibuprofen	55	1 - 7	No interference
Lactose	1000	< 0.5	No interference
Lidocaine	6	0.15 - 0.6	No interference
Magnesium	5 mM	1.1 mM	No interference
Maltose	1000	N/A	No interference
Mannitol	1000	0.0128	No interference
Metaproterenol	1.81	0.00022 - 0.00130	No interference
Metoprolol	0.3	0.005 - 0.027	No interference
Naproxen	100	3-12	No interference
Nifedipine	0.17	0.017	No interference
Nortriptylline	0.15	0.005 - 0.015	No interference
Penicillin	12	1.2	No interference
pH	6.7 - 9.8	7.35 - 7.45	No interference
Phenytoin	10	1 - 2	No interference
Piroxicam	5	0.3-0.5	No interference
Potassium	10 mM	5.9 mM	No interference
Sodium	200 mM	135 - 145 mM	No interference
Sorbitol	1000	0.044	No interference
Sulfamethoxazole	120	5-12	No interference
Sulfate	5 mM	1 mM	No interference
Terfenadine	0.45	0.00015 - 0.00045	No interference
Tetracycline	4	0.4	No interference
Theophylline	25	1.0 - 2.0	No interference
Tolbutamide	64	4.32 - 24	No interference
Total Protein (gamma-Globulin)	12000	6000 - 8000	No interference
Triglycerides	2000	30-300	No interference
Urea	600	38	No interference
Vancomycin	25	0.025	No interference
Verapamil	0.45	0.014 - 0.045	No interference
Vitamin E	20	0.5 - 2.0	No interference
Warfarin	2	0.1 - 1.0	No interference

FORA GD20

Blood Glucose Monitoring System



For additional product information,
contact Links Medical Products
at (888) 425-1149.



Links Medical Products Inc.[®]
www.linksmed.com