

JXB-183

Non-Contact Infrared Thermometer

Device Type	Non-contact infrared thermometer
Measurement Distance	1.2" to 2" (3 to 5 cm) from forehead
Measuring Range	Body temperature: 89.6°F to 109.4°F (32°C to 43°C) Surface temperature: 32°F to 140°F (0°C to 60°C) Room temperature: 32°F to 104°F (0°C to 40°C)
Display Type / Resolution	Large, 3-color LCD with highly accurate 1/10 degree reading
Accuracy	89.6°F to 94.9°F=±0.6°F (32°C to 34.9°C=±0.3°C) 95°F to 107.6°F=±0.4°F (35°C to 42°C=±0.2°C) 107.7°F to 109.4°F=±0.6°F (42.1°C to 43°C=±0.3°C)
Normal Using Condition (Ambient Temperature)	50°F to 104°F (10°C to 40°C); Relative humidity ≤ 85%
Consumption	≤300 mW
Calibration	Preset at factory. Device re-calibrates each time device is turned ON.
Memory	Stores up to 32 readings
Automatic Power Off	30 seconds after last measurement
Batteries	2 AAA batteries (not included)
Unit Weight	67g without batteries
Warranty	1 Year limited warranty (manufacturer defect only)



For more information about the JXB-183 Non-Contact Infrared Thermometer, please contact your local distributor or Links Medical Products customer service at (888) 425-1149.

LMP Item# LMP005

15979 North 76th Street, Unit D
Scottsdale, AZ 85260 USA
www.linksmed.com

LMP
Links Medical Products Inc.®