

JXB-183

Non-Contact Infrared Thermometer

HELPS PREVENT
THE SPREAD OF
INFECTIOUS DISEASE

Provides quick & accurate measurement
of core body temperature

- No contact with resident's skin reduces the risk of cross-contamination
- FDA approved, clinical-grade thermometer for multi-patient and individual use
- No probe covers required - significant annual savings
- Ideal for sleeping or combative residents
- Ideal for F-Tag 257 compliance (comfortable & safe temperature levels)

ACCURATE
and precise

INSTANT
1 second
measurement

AUDIBLE
option



FEATURES

- Non-contact reading in one second from 1.2 - 2 inches distance from forehead
- LCD digital display screen with high temperature alarm
- Stores the previous 32 temperature measurements
- Temperature can be displayed in either °F or °C
- Automatic power-off (30 seconds) to conserve energy
- Modes for body, surface and room (ambient) temperature
- Visual/audible alarm if temperature exceeds 100.4°F
- 1 Year limited warranty - manufacturer defect only


FDA APPROVED

LMP Item# LMP005

15979 North 76th Street, Unit D
Scottsdale, AZ 85260 USA
Toll free: (888) 425-1149
Fax: (949) 753-7412
www.linksmed.com


Links Medical Products Inc.®

VALIDATED ACCURACY

The JXB-183 non-contact thermometer has a validated accuracy of +/- 0.2 °C (0.4 °F) over the clinical range of valid temperature.*

* Standard for ISO 13485