

Spec Sheet

LINKTEMP[®]

Non-Contact Infrared Thermometer

Device Type	Non-contact infrared thermometer
Measurement Distance	1.2" to 2" (3 to 5 cm) from forehead
Measuring Range	Body temperature: 89.6°F to 109°F (32°C to 42.9°C) Surface temperature: 32°F to 140°F (0°C to 60°C) Ambient temperature: 32°F to 104°F (0°C to 40°C)
Displayed Temperature Range	89.6°F to 109.2°F (32°C to 42.9°C)
Display Type	Large, 2-color LCD with highly accurate 1/100 degree reading
Display Resolution	0.1°F (0.1°C)
Accuracy	89.6°F to 94.8°F=±0.6°F (32°C to 34.9°C=±0.3°C) 95°F to 107.6°F=±0.4°F (35°C to 42°C=±0.2°C) 107.8°F to 109.2°F=±0.6°F (42.1°C to 42.9°C=±0.3°C)
Operating Temperature	50°F to 104°F (10°C to 40°C)
Operating Humidity	Up to 85% non condensing
Storage Temperature	-4°F to 131°F (-20°C to 55°C)
Storage Humidity	Up to 95% non condensing
Pressure Altitude	700 hPa - 1060 hPa
Consumption	≤450 mW
Calibration	Preset at factory. Device re-calibrates each time device is turned ON
Memory	Stores up to 32 readings
Automatic Power Off	30 seconds after last measurement
Batteries	2 AA batteries
Unit Weight	3.48 oz without batteries
Unit Size	6.81" x 2.05" x 1.77" (173 x 52 x 45 mm) in L x W x H
Warranty	18 month limited warranty



For more information on the **LinkTemp™** Non-Contact Infrared Thermometer, please contact your local distributor or Links Medical Products customer service at (888) 425-1149.



Links Medical Products Inc.[®]

15979 North 76th Street, Unit D
Scottsdale, AZ 85260 USA
www.linksmed.com