

# LINKTEMP<sup>®</sup>

## Non-Contact Infrared Thermometer

**Audible  
response in  
both English &  
Spanish**

**Provides quick & accurate measurement  
of core body temperature**

- ◇ No probe covers required - significant annual savings
- ◇ No contact with resident's skin reduces the risk of cross-contamination
- ◇ Ideal for sleeping or combative residents
- ◇ FDA cleared, clinical-grade thermometer for multi-patient and individual use
- ◇ Ideal for F-Tag 257 compliance (comfortable & safe temperature levels)

## FEATURES

- Non-contact reading in one second from 1.2 - 2 inches distance from forehead
- Two-color LCD with highly accurate 1/100 degree reading
- Stores the previous 32 temperature measurements
- Temperature can be displayed in either °F or °C
- Automatic power-off (30 seconds) to conserve energy
- Modes for body, surface and room (ambient) temperature
- Visual/audible alarm if temperature exceeds 100.4 °F
- Performs up to 100,000 readings
- Long battery life

**ACCURATE**  
and precise

**INSTANT**  
1 second  
measurement

**AUDIBLE**  
bilingual voice  
option



LMP Item # LMP001  
2 AA batteries included  
18 month limited warranty

## VALIDATED ACCURACY

The LinkTemp non-contact thermometer has a validated accuracy of +/- 0.2 °C (0.4 °F) over the clinical range of valid temperature.\*

\* Standard for ISO 13485

15979 North 76th Street, Unit D  
Scottsdale, AZ 85260 USA  
Toll free: (888) 425-1149  
Fax: (949) 753-7412  
www.linksmed.com

**lmp**  
Links Medical Products Inc.®