

SILENT KNIGHT®

Advanced Pill Crushing Device



Protected by US patents: 7,051, 963 B1, D497, 543 S, 5, 915,637, 9,125,796.

ACCLAIMED BENEFITS OF THE NEW SILENT KNIGHT®

- SILENT** Won't disturb patients or staff.
- ERGONOMIC HANDLE** Fits all hand sizes.
- EFFECTIVE** Powerful dual action crushes pills to powder in a few seconds.
- PATIENT SAFETY** Eliminates cross contamination.
- USER SAFETY** Maximum results with minimal effort. Contour cover encloses internal moving parts.

EXCLUSIVE PILL CRUSHER POUCH

The Silent Knight® does not function alone. Its companion item the patented Silent Knight® Pill Crusher Pouch is essential to guarantee maximum performance while reducing cost, waste, malfunction and lost time.

The pouch's distinctive curved bottom structure assures full recovery of crushed medication. The thick, durable construction effectively reduces pharmaceutical cross contamination and eliminates loss of medication.



Note: No other pill crusher pouch will effectively and safely function in the Silent Knight®. Use of any other pill pouch will compromise product performance and outcome.

ADVANCED ENGINEERING & DESIGN

Advanced engineering and design components have been applied to the original Silent Knight's basic operating and crushing features to create the New Silent Knight.

The patented dual action "crushing and grinding" movements of the New Silent Knight result in silent, smooth, effortless and effective crushing of pills to fine powder.

The crusher handle is designed for all hand sizes to fit comfortably. User safety is enhanced with the specially designed internal handle shock absorber. And the surgical quality stainless steel crushing components, along with a solid aluminum base and handle, make the New Silent Knight a powerful, highly effective and durable piece of medical equipment.



Ribbed crusher plate action

ORDERING INFORMATION:

Silent Knight Pill Crusher:
Item Number SK0500

Silent Knight Pill Crusher Pouch:
Item Number PC1000

WARRANTY:

Limited Lifetime Warranty

Contact your preferred distributor to purchase the new Silent Knight Pill Crusher.

You can also contact the Links Medical customer service department at (888) 425-1149 for a list of distributors in your area that stock the new Silent Knight Pill Crusher.